

UROMED TBK-DUFOUR

REF 1870/1871

Dufourspitze

Zweilumig, 2 Augen, Länge 43 cm

Tamponade-Ballon-Katheter zur postoperativen Blasendrainage

- lange, knickstabile Spitze
- hohe Drainageleistung
- röntgenologische Darstellbarkeit der Spitze
- große Augen
- reines Silikon
- fest

Tamponade Balloon Catheter for postoperative bladder drainage

- 2-way, 2 eyes, length 43 cm
- long, nonkinking tip
- high drainage results
- radiopaque tip
- large eyes
- pure Silicone
- solid



REF REF	Größen Sizes	Ballon Balloon	PZN PZN
1870 - 18	Ch. 18	40 ml	04672707
1870 - 20	Ch. 20	50 ml	04672742
1870 - 22	Ch. 22	50 ml	04672771
1870 - 24	Ch. 24	50 ml	04672802

Hilfsmittelkatalog: 15.25.15.6033

REF REF	Größen Sizes	Ballon Balloon	PZN PZN
1871 - 20	Ch. 20	60 ml	04672825
1871 - 22	Ch. 22	60 ml	04672854

Hilfsmittelkatalog: 15.25.15.6033

Verpackung: doppelt, steril
Original-Karton: 10 Stück

Packing: double, sterile
Original carton: 10 pieces

UROMED
PRODUKTE FÜR DIE UROLOGIE