

# UROMED URETHROTOMIE BALLON-KATHETER

## REF 3200

### Zylindrische Spitze

Dreilumig, 2 Augen, 3 Löcher  
Länge 43 cm

Durch die Perforation kann nach Urethrotomia interna sowohl ein Sekretabfluss als auch eine Medikamenteninstillation erfolgen.

*The perforation allows the drainage of secretion as well as the instillation of drugs after urethrotomy interna.*

### Ballon-Katheter zur transurethralen Medikamentenapplikation / Blasen- drainage

- Perforation im Bereich der bulbären Urethra
- großlumig
- röntgenologische Darstellbarkeit der Spitze
- reines Silikon
- weich

### *Ballon Catheter for the transurethral application of drugs / bladder drainage*

- *3-way, 2 eyes, 3 holes, length 43 cm*
- *perforation in the area of the bulbar urethra*
- *wide lumen*
- *radiopaque tip*
- *pure Silicone*
- *soft*



REF REF	Größen Sizes	Ballon Balloon
3200 - 16	Ch. 16	5 – 10 ml
3200 - 18	Ch. 18	5 – 10 ml
3200 - 20	Ch. 20	5 – 10 ml
3200 - 22	Ch. 22	5 – 10 ml
3200 - 24	Ch. 24	5 – 10 ml

**Verpackung:** doppelt, steril  
**Original-Karton:** 10 Stück

**Packing:** double, sterile  
**Original carton:** 10 pieces