

UROMED SUPRAPUBISCHER BALLON-KATHETER

REF 2956/2966


Kurze, zentral offene Spitze
2 Augen

Ballon-Katheter zur suprapubischen Blasendrainage

- kurze, atraumatische Spitze mit abgerundeter Kante
- flache Ballonform
- großlumig
- Soft-Ballon
- reines Silikon
- weich

Ballon Catheter for suprapubic drainage of the urinary bladder

- short, atraumatic and central open tip with rounded edge, 2 eyes
- flat shape of balloon
- wide lumen
- soft balloon
- pure silicone
- soft



REF 2956
Länge 43 cm
length 43 cm

REF 2966
Länge 33 cm
length 33 cm

REF REF	Größen Sizes	Ballon Balloon	PZN PZN
2956-10	Ch. 10	2 ml	00051724
2956-12	Ch. 12	5 ml	00051730
2956-14	Ch. 14	5 ml	00051747
2956-16	Ch. 16	5 ml	00051753
2956-18	Ch. 18	5 ml	00051776
2956-20	Ch. 20	5 ml	00051807
2956-22	Ch. 22	5 ml	02784388
2956-24	Ch. 24	5 ml	02784394

REF REF	Größen Sizes	Ballon Balloon	PZN PZN
2966-10	Ch. 10	2 ml	07654443
2966-12	Ch. 12	5 ml	04741957
2966-14	Ch. 14	5 ml	03281465
2966-16	Ch. 16	5 ml	04741963
2966-18	Ch. 18	5 ml	04741986
2966-20	Ch. 20	5 ml	04741992
2966-22	Ch. 22	5 ml	02784402
2966-24	Ch. 24	5 ml	02784431

Verpackung: doppelt, steril
Original-Karton: 10 Stück

Packing: double, sterile
Original carton: 10 pieces

UROMED
PRODUKTE FÜR DIE UROLOGIE