

UROMED »HYDROPUR®« STEUERBAR

UROMED »HYDROPUR®« MANOEUVRABLE

REF 4485 – 4488

**Zylindrische Spitze
mit Perforation im Ureterteil**
*Closed tip with perforation
in the ureteral part*

**Katheter zur inneren
Schienung des Ureters**

Für die ureterorenoskopische
Anwendung ist die Ureter-
Dauerschleife in Ch. 4,8 vorgesehen.

- aromatisches Polyurethan
- mit vergüteter Oberfläche
- geringe Inkrustationsneigung
- zur Langzeitanwendung geeignet
- gutes Vorschiebehaviorn durch mittlere Materialkonsistenz
- konische, dilatierende Spitze
- kleine, kreisrunde Pigtails
- atraumatische Augen
- vormontiertes Set zur einfachen Handhabung
- röntgenologisch darstellbar
- graduiert
- steuerbar

**Catheter for the internal stenting
of the ureter**

*The ureteral stent in Ch. 4,8 is destined
for ureterorenoscopic application.*

- aromatic polyurethane
- with refined surface
- low incrustation rate
- suitable for long-term application
- good advancement control because of medium material consistency
- conical, dilating tip
- small-sized, circular pigtails
- atraumatic eyes
- pre-assembled set for easy handling
- overall radiopaque
- graduated
- manoeuvrable



Ureter-Dauerschleife
beidseitiger Pigtail

*ureteral stent, double pigtail
cylindrical tip, with perforation
in the ureteral part*

Ch. 6-Ch. 8
Pushersystem
Länge 45 cm

*Ch. 6-Ch. 8
pusher system
length 45 cm*

Ch. 4,8
Pushersystem
Länge 70 cm

*Ch. 4,8
pusher system
length 70 cm*

Spiralmandrin
PTFE-beschichtet
guidewire
PTFE coated

1 Fixierklemme bei REF 4485
1 fixing clamp for REF 4485

REF REF	Dauerschleife Stent	Mandrin Guidewire	PZN PZN
4485 - 24	Ch. 4,8 / 24 cm	0,028" / 150 cm	11371739
4485 - 26	Ch. 4,8 / 26 cm	0,028" / 150 cm	11371745
4485 - 28	Ch. 4,8 / 28 cm	0,028" / 150 cm	11371751
4485 - 30	Ch. 4,8 / 30 cm	0,028" / 150 cm	11371768
4486 - 24	Ch. 6,0 / 24 cm	0,035" / 100 cm	11371774
4486 - 26	Ch. 6,0 / 26 cm	0,035" / 100 cm	11371780
4486 - 28	Ch. 6,0 / 28 cm	0,035" / 100 cm	11371797
4486 - 30	Ch. 6,0 / 30 cm	0,035" / 100 cm	11371805

REF REF	Dauerschleife Stent	Mandrin Guidewire	PZN PZN
4487 - 24	Ch. 7,0 / 24 cm	0,035" / 100 cm	11371811
4487 - 26	Ch. 7,0 / 26 cm	0,035" / 100 cm	11371828
4487 - 28	Ch. 7,0 / 28 cm	0,035" / 100 cm	11371834
4487 - 30	Ch. 7,0 / 30 cm	0,035" / 100 cm	11371840
4488 - 26	Ch. 8,0 / 26 cm	0,038" / 100 cm	11371857
4488 - 28	Ch. 8,0 / 28 cm	0,038" / 100 cm	11371863
4488 - 30	Ch. 8,0 / 30 cm	0,038" / 100 cm	11371886

Verpackung: doppelt, steril
Original-Karton: 10 Stück

Packing: double, sterile
Original carton: 10 pieces

UROMED
PRODUKTE FÜR DIE UROLOGIE