

UROMED Führungsdraht

UROMED Guide Wire

REF 4271 - 4295

Titandraht,
hydrophil beschichtet
titanium wire,
hydrophilic coated

Als „Pfadfinder“ hat sich dieses Instrument bei der Überwindung schwierigster Passagen bewährt.



This instrument has proven its worth as "pathfinder" even in most difficult passages.

Flexibles Führungsinstrument für urologische Behandlungen

- atraumatische Spitze
- hydrophile Beschichtung
- extrem gleitfähig
- dreh- und knickstabil
- flexibel und hochelastisch
- minimiert das Perforationsrisiko
- vor Gebrauch Aktivierung der hydrophilen Beschichtung im Dispenser
- gerade und abgewinkelte Ausführung
- auch in extra steifer Ausführung erhältlich
- rund verpackt

Flexible guidance instrument for urological treatments

- atraumatic tip
- hydrophilic coating
- extremely high sliding ability
- nontwisting and nonkinking
- flexible and highly elastical
- minimizes the risk of perforation
- activating the hydrophilic coating in the dispenser before use
- straight and angled tip
- available in extra stiff
- round packaging



REF 4271 - 4275

abgewinkelte Spitze
angled tip

REF 4290

abgewinkelte Spitze, extra steif
angled tip, extra stiff

REF 4281 - 4285

gerade Spitze
straight tip

REF 4295

gerade Spitze, extra steif
straight tip, extra stiff

REF <i>REF</i>	Länge <i>Length</i>	Spitzen <i>Tips</i>	Durchmesser <i>Diameters</i>	Steifigkeit <i>Stiffness</i>
4271	150 cm	abgewinkelt/ <i>angled</i>	0,025"/0,64 mm	standard/ <i>regular</i>
4273	150 cm	abgewinkelt/ <i>angled</i>	0,032"/0,82 mm	standard/ <i>regular</i>
4275	150 cm	abgewinkelt/ <i>angled</i>	0,035"/0,89 mm	standard/ <i>regular</i>
4281	150 cm	gerade/ <i>straight</i>	0,025"/0,64 mm	standard/ <i>regular</i>
4283	150 cm	gerade/ <i>straight</i>	0,032"/0,82 mm	standard/ <i>regular</i>
4285	150 cm	gerade/ <i>straight</i>	0,035"/0,89 mm	standard/ <i>regular</i>
4290	150 cm	abgewinkelt/ <i>angled</i>	0,035"/0,89 mm	extra steif/ <i>extra stiff</i>
4295	150 cm	gerade/ <i>straight</i>	0,035"/0,89 mm	extra steif/ <i>extra stiff</i>

UROMED
PRODUKTE FÜR DIE UROLOGIE

Verpackung: steril
Originalkarton: 5 Stück

Packaging: sterile
Original carton: 5 pieces