

# UROMED UROject Injektionskanüle

## UROMED UROject Injection Cannula

**REF 6240**

**Stichlänge 5 mm**  
**Schaftlänge 650 mm**


**Needle length 5 mm**  
**Length of shaft 650 mm**

### Flexible Injektionskanüle zur Botox®-Injektion in die Blasenwand

- individuell für starres und für flexibles Zystoskop geeignet
- gute Führbarkeit durch distale Materialverstärkung
- exzellente Dosierbarkeit von Botox® durch Leervolumenangabe und sehr feine Nadel
- gute Haptik sowie sicheres Aufsetzen und Entfernen der Spritze durch Luer-Lock-Verbindung
- Nadelschutz für das sichere Einführen in das Instrument
- gleichbleibend scharfer Nadelanschiff gewährleistet eine präzise und sichere Injektion

### Flexible injection cannula for Botox® injection into the bladder wall

- suitable for rigid as well as for flexible cystoscopes
- reinforcement of distal material made for good guidance
- excellent Botox® dosing because of empty volume specification and very fine needle
- good haptics, safe attachment and removal of the syringe because of luer-lock connection
- needle protection for safe insertion into the instrument
- consistently sharp bevelled needle guarantees accurate and safe injection



REF	Ch./Fr.	Ø Schaft	Schaftlänge	Ø Kanüle	Stichlänge	PZN
REF	Ch./Fr.	Ø Shaft	Length of shaft	Ø Cannula	Needle length	CPN
6240	05	1,7 mm	650 mm	27 G/0,4 mm	5 mm	12364581

**Verpackung:** steril  
**Originalkarton:** 10 Stück

**Packaging:** sterile  
**Original carton:** 10 pieces