

# UROMED CYSTOMETRIE-KATHETER **PURgreen**<sup>®</sup>

## UROMED CYSTOMETRY CATHETER **PURgreen**<sup>®</sup>

### REF 2934

male, zylindrische Spitze  
dreilumig, Länge 40 cm


*male, cylindrical tip  
3-way, length 40 cm*

### Katheter zur transurethralen Druckmessung in der Harnblase/Urethra

- für Geräte mit Rollenpumpe
- je 1 Auge pro Messkanal
- Augenzahl 1+2 (1 Auge an Spitze)
- großlumiger Zulauf mit 1 Auge
- präzise Messergebnisse
- Luer-Lock-Anschlüsse mit Kappe
- röntgenologisch darstellbar
- graduiert
- reines Polyurethan
- halbfestes, flexibles Material

### Catheter for transurethral measurement of pressure in the urinary bladder/urethra

- for instruments with roller pump
- 1 eye each per measuring channel
- 1+2 eyes (1 eye at tip)
- wide lumen inflow with 1 eye
- precise measuring results
- Luer-Lock adapters with cap
- overall radiopaque
- graduated
- pure polyurethane
- semi-solid, flexible material



REF	Größe
REF	Size
2934.09	Ch. 9/3,0 mm

**Verpackung:** doppelt, steril  
**Originalkarton:** 10 Stück

**Packaging:** double, sterile  
**Original carton:** 10 pieces