

UROMED CYSTOMETRIE-KATHETER

UROMED CYSTOMETRY CATHETER

REF 2923

female, zylindrische Spitze
zweilumig, Länge 20 cm

*female, cylindrical tip
2-way, length 20 cm*

Das UROMED Cystometrie-Programm PURgreen aus reinem Polyurethan zeichnet sich durch seine hochwertige Qualität und präzise Messergebnisse aus. Die zweilumigen Modelle »female« und »male« weisen unterschiedliche, den Erfordernissen angepasste Positionierungen der Augen auf. Produkte und Verpackung sind PVC- und latexfrei.


The high quality and precise measuring results distinguish the UROMED cystometry programme PURgreen made of pure polyurethane. The 2-way versions »female« and »male« show different positionings of the eyes which are adapted to the requirements. Products and packaging are free of PVC and latex.

Katheter zur transurethralen Druckmessung in der Harnblase

- Messkanal mit 2 Augen
- großlumiger Zulauf mit 2 Augen
- präzise Messergebnisse
- Luer-Lock-Anschlüsse mit Kappe
- röntgenologisch darstellbar
- graduert
- reines Polyurethan
- halbfestes, flexibles Material

Catheter for transurethral measurement of pressure in the urinary bladder

- measuring channel with 2 eyes
- wide lumen inflow with 2 eyes
- precise measuring results
- Luer-Lock adapters with cap
- overall radiopaque
- graduated
- pure polyurethane
- semi-solid, flexible material



REF	Größen
REF	Sizes
2923.07	Ch. 7 / 2,3 mm
2923.09	Ch. 9 / 3,0 mm

Verpackung: doppelt, steril
Originalkarton: 10 Stück

Packaging: double, sterile
Original carton: 10 pieces